

BS-821 HOSPITAL BED



High Quality

1. Advanced medical motor system, including 2 pcs electric motors, 1 pcs control box, 1 pcs Cable remote control
2. 1 couple of rail sides, made of aluminium alloy, good for anti-rust, High Hardness with long time life 4 PCS Infusion Urine hooks, and 1 pcs telescopic SS IV pole.
3. 125mm dia. Noiseless castors with brake, stable standing and smooth moving for common ground Single touch style handle control, convenient to control to bed for nurse

01

We strictly control the details of each product. Never show any flaws.



02

Design simple, convenient, generous, With a breath of life.



03

Intelligent adjustment height, a variety of stalls, a variety of angles, a key operation.



04

Let you sleep comfortable, peace of mind, good mood, quick recovery .

